

## Universal IP/Co-IP Toolkit (Agarose)

Cat #: KTD105-EN

Size: 20 T

	<b>Universal IP/Co-IP Toolkit (Agarose)</b>		
	<b>Cat #:</b> KTD105-EN		<b>Lot #:</b> Refer to product label
	<b>Applicable samples:</b> Animal and Plant Tissues, Cells		
	<b>Storage:</b> Store according to the recommended storage conditions of each component, stable for 12 months		

## Assay Principle

Immunoprecipitation (Immunoprecipitation, IP) is a small affinity purification method using specific antibodies fixed on solid-phase supports such as magnetic beads or agarose. As an important part of many proteomics-related studies, IP can be used to detect the existence of proteins, relative abundance, up-and-down expression of proteins, stability and interaction of proteins and so on. But they are tedious and time-consuming, with a wide variety of reagents, and poor stability and reliability of experimental results. Based on the pain points of the traditional IP/Co-IP experiment, Abbkine developed a Universal IP/Co-IP Toolkit (Agarose); After designed and researched, the toolkit contains: ① The optimized natural and denatured lysate can meet the needs of IP/Co-IP or WB at the same time; ② Compared with Protein A or Protein G Agarose, Protein A/G Agarose has a wider binding range, more efficient antibody binding capacity and lower protein non-specific adsorption rate, making the operation more time-saving, efficient and saving antibody consumption. ③ The ready-to-use IP negative control can eliminate the non-specific binding of IgG itself to the target protein or other specific biomolecules, ensure the specificity of IP antibodies, and make the experimental data more stable and reliable. ④ The unique IPkine™ second antibody (HRP labeled goat anti-mouse IgG light chain secondary antibody and HRP labeled mouse anti-rabbit IgG light chain secondary antibody) can perfectly eliminate heavy chain interference. Universal IP/Co-IP Toolkit (Agarose), as a standardized product of IP/Co-IP, can meet the IP/Co-IP needs of most users, and it's the best choice for IP/Co-IP detection experiments. In addition, the corresponding products (A25012, A25022) in the Toolkit can also be ordered separately.

## Materials Supplied and Storage Conditions

Toolkit components	Kit components	Size	Storage conditions
Universal IP/Co-IP Toolkit (Agarose) (Part 1 of 2)	Non-Denaturing Lysis Buffer	25 mL	4°C
	Denaturing Lysis Buffer	25 mL	4°C
	10×Wash Buffer	20 mL	RT
	Protein A/G Agarose	0.45 mL	4°C
	Elution Buffer	2 mL	4°C

	Neutralization Buffer	0.2 mL	4°C
Universal IP/Co-IP Toolkit (Agarose) (Part 2 of 2)	100×Proteinase Inhibitor Cocktail	0.2 mL	-20°C
	Mouse IgG (1mg/mL)	30 μL	-20°C
	Rabbit IgG (1mg/mL)	30 μL	-20°C
	IPKine™ HRP, Goat Anti-Mouse IgG LCS	30 μL	-20°C
	IPKine™ HRP, Mouse Anti-Rabbit IgG LCS	30 μL	-20°C

## Materials Required but Not Supplied

- Vertical rotating mixer
- Freezing Centrifuge
- Precision Pipettes, Disposable Pipette Tips
- Deionized Water
- PBS Buffer
- Dounce homogenizer (for tissue samples)
- SDS-PAGE Loading Buffer

## Reagent Preparation

**Non-Denaturing Lysis Buffer:** Native protein lysis buffer, extract protein for IP samples. Ready to use as supplied. Place it on ice for use. Store at 4°C.

**Denaturing Lysis Buffer:** Denatured protein lysis buffer, extract protein for WB samples. Ready to use as supplied. Place it on ice for use. Store at 4°C.

**1×Wash Buffer:** Add Deionized Water to the 10× Buffer and dilute the 10×Wash Buffer to 1×Wash Buffer before use. Store at 4°C.

**Protein A/G Agarose:** Ready to use as supplied. Store at 4°C, Avoid frozen.收

**Note:** 10×Wash Buffer stored at 4°C will produce precipitation, and can be preheated in 37°C water bath for 10 min before use to dissolve precipitation.

**Elution Buffer:** Ready to use as supplied. Elute Non-Denaturing protein. Store at 4°C.

**Neutralization Buffer:** Ready to use as supplied. Neutralize Non-Denaturing eluted protein. Store at 4°C.

**100×Proteinase Inhibitor Cocktail:** Add 1×Proteinase Inhibitor Cocktail to Lysis Buffer before use. Store at -20°C.

**Mouse IgG:** Ready to use as supplied. Used as a negative control for Mouse IP antibody. Store at -20°C.

**Rabbit IgG:** Ready to use as supplied. Used as a negative control for Rabbit IP antibody. Store at -20°C.

**IPKine™ HRP, Goat Anti-Mouse IgG LCS:** HRP conjugated Goat Anti-Mouse IgG light chain specific secondary antibody. The recommended dilution ratio of antibody for WB is 1:1000-1:10,000 (optimal 1:2,000). Store at -20°C.

**IPKine™ HRP, Mouse Anti-Rabbit IgG LCS:** HRP conjugated Mouse Anti-Rabbit IgG light chain specific secondary antibody. The recommended dilution ratio of antibody for WB is 1:1000-1:10,000 (optimal 1:2,000). Store at -20°C.

## Assay Procedure

### I Immunoprecipitation

#### A. Preparation of protein samples

##### 1. Extract protein for Cell Samples:

(1) Collect cells (Adherent cells: 80% to 90% of monolayer cells were grown in a 10 cm cell culture dish. Remove the medium and wash with PBS once; Suspended cells: Collect  $5 \times 10^6$  cells by centrifugation and washed by PBS once)

(2) Add 0.5-1 mL ice-cold Non-Denaturing Lysis Buffer to cells (Add 1×Proteinase Inhibitor Cocktail to Non-Denaturing Lysis

Buffer), lytic cells at 4°C for 5 min. During the process, the pipette is used to blow the mixture repeatedly, transfer cell suspension to a new tube.

(3) Centrifuge at 12,000 rpm for 10 min at 4°C, collect supernatant, and detect protein concentration by BCA Assay (It is recommended to use Abbkine Cat #: KTD3001 Protein Quantification Kit (BCA Assay)).

## 2. Extract protein for Tissue Samples:

(1) Plant or Animal Tissue Samples: Weigh 0.1 g of tissue and add 1 mL Non-Denaturing Lysis Buffer (Add 1×Proteinase Inhibitor Cocktail to Non-Denaturing Lysis Buffer), Homogenize tissue with Liquid nitrogen or Dounce homogenizer. (If the protein concentration is low, reduce the volume of Non-Denaturing Lysis Buffer)

(2) Transfer the homogenate to a new tube, lytic samples at 4°C for 5 min.

(3) Centrifuge at 12,000 rpm for 10 min at 4°C, collect supernatant, and detect protein concentration by BCA Assay. (It is recommended to use Abbkine Cat #: KTD3001 Protein Quantification Kit (BCA Assay))

**Note: The optimal concentration of total protein was in the range of 0.5-1 µg/µL. Usually, the total protein concentration needs to be pre-adjusted according to the different expression levels of the target protein. Co-Immunoprecipitation (Co-IP) usually requires the use of fresh samples that have not been frozen. Although frozen protein samples can be used for immunoprecipitation, fresh samples are preferred.**

### B. Remove non-specifically binding (Optional):

(1) Add 20 µL (40 µL suspension) Protein A/G Agarose to a 1.5 mL centrifuge tube, Centrifuge at 800 rpm for 2 min at 4°C, remove the supernatant.

**Note: Protein A/G Agarose must be fully suspended before use, that is, fully invert the mixture several times to mix well. If Protein A/G Agarose agglomerated which is caused by the reduction of protective solution (0.1 M, pH 7.4 PBS containing 20% ethanol), it is recommended to add protective solution, and keep the same volume as Agarose. For per mL Protein A/G Agarose contains 0.5 mL Agarose and 0.5 mL protective solution, each sample need add 40 µL Protein A/G Agarose suspension.**

(2) Add 1mL 1×Wash Buffer, and re-suspend Protein A/G Agarose, Centrifuge at 800 rpm for 2 min at 4°C, remove the supernatant, repeat 3 times.

(3) Add 0.2-1 mL(0.2-1 mg) protein sample, shake and incubate at 4°C for 30 min.(It is recommended to use vertical rotating mixer with Low-speed rotation)

(4) Centrifuge at 800 rpm for 2 min at 4°C, use supernatant for Immunoprecipitation.

**Note: During centrifugation, the centrifugal speed should not be too high, which will damage Protein A/G Agarose. The centrifugal speed should be 800-2,000 rpm.**

### C. Immunoprecipitation

(1) Add 20 µL (40 µL suspension) Protein A/G Agarose to a 1.5 mL centrifuge tube, centrifuge at 800 rpm for 2 min at 4°C, remove the supernatant.

(2) Add 1mL 1×Wash Buffer, and re-suspend Protein A/G Agarose, Centrifuge at 800 rpm for 2 min at 4°C, remove the supernatant, repeat 3 times.

(3) Add 0.2-2 µg antibodies, and re-suspend Protein A/G Agarose, shake and incubate at room temperature for 30 min. Centrifuge at 800 rpm for 2 min, collect the supernatant, use the precipitation for Step (4). **Optional, IgG negative control: Add 0.2-2 µL Mouse IgG or Rabbit IgG, and re-suspend Protein A/G Agarose, shake and incubate at room temperature for 30 min. Centrifuge at 800 rpm for 2 min, collect the supernatant, use the precipitation for Step (4). This step could exclude the non-specific binding of IgG to the target protein or other specific biological molecules.**

**Note: Protein A/G Agarose needs fluid to keep it flowing while use vertical rotating mixer, if it does not flow, add a certain volume of 1×Wash Buffer.**

(4) Add 1 mL 1×Wash Buffer, and re-suspend Protein A/G Agarose, centrifuge at 800 rpm for 2 min at 4°C, remove the supernatant, repeat 3 times.

(5) Add 0.2-1 mL(0.2-1 mg) protein sample or supernatant (supernatant which remove non-specifically binding), shake and incubate at room temperature for 1 h or at 4°C overnight.

- (6) Centrifuge at 800 rpm for 2 min at 4°C, collect the supernatant, use the precipitation for Step (7).
- (7) Add 1 mL 1×Wash Buffer, and re-suspend Protein A/G Agarose, Centrifuge at 800 rpm for 2 min at 4°C, remove the supernatant, repeat 3 times.
- (8) Elution of the Immune Complex.
- a) Denatured elution: This method is suitable for SDS-PAGE and Western Blotting analysis of elution samples. Add 20-50 µL 1×SDS-PAGE Loading Buffer to the tube and mix well, incubate at 100°C for 5 min, then Centrifuge at 800 rpm for 1 min, and collect the supernatant to a new tube for SDS-PAGE and Western Blotting analysis.
- b) Non-Denatured elution: The eluted samples retained their original biological activity and could be used for subsequent functional analysis. Add 20-50 µL Elution Buffer to the tube and mix well, and incubate at room temperature for 10 min, then centrifuge at 800 rpm for 2 min at 4°C. Take the supernatant to collect the elution component, which is the target antigen. Transfer the supernatant to a new tube, and immediately add 1/10 Elution Buffer volume Neutralization Buffer to adjust the pH to 7.0-8.0, use the elution sample for subsequent functional analysis.

## II Co-Immunoprecipitation

Refer to the steps of Immunoprecipitation.

**Strawberry moment:** The protein with high purity can be obtained with Universal IP/Co-IP Toolkit conveniently, and it also can be obtained with different tag mouse monoclonal antibody (Cat #: A02010AGB, A02040AGB, A02050AGB) by Immunoprecipitation or Co-Immunoprecipitation. Scan the QR code on the right and follow the Abbkine official account to learn more about Abbkine products.



## Recommended Products

Catalog No.	Product Name	Recommended reason
KTP2070	PurKine™ Antibody Purification Kit (Protein A/G)	High throughput, high performance, high stability
A25012	IPKine™ HRP, Goat Anti-Mouse IgG LCS	Specifically identify the light chain of IgG and avoid the interference of heavy chain
A25022	IPKine™ HRP, Mouse Anti-Rabbit IgG LCS	Specifically identify the light chain of IgG and avoid the interference of heavy chain
A25112	IPKine™ HRP, Goat Anti-Mouse IgG HCS	Specifically identify the heavy chain of IgG and avoid the interference of light chain
A25122	IPKine™ HRP, Mouse Anti-Rabbit IgG HCS	Specifically identify the heavy chain of IgG and avoid the interference of light chain
ABT2013	Agarose Conjugated Anti-DDDDK Tag Mouse Monoclonal Antibody (1B10)	Featured DDDDK Tag Mouse Monoclonal Antibody, designed for immunoprecipitation (IP) of tagged protein
ABT2043	Agarose Conjugated Anti-HA Tag Mouse Monoclonal Antibody (4F6)	Featured HA Tag Mouse Monoclonal Antibody, designed for immunoprecipitation (IP) of tagged protein
ABT2063	Agarose Conjugated Anti-Myc Tag Mouse Monoclonal Antibody (2D5)	Featured Myc Tag Mouse Monoclonal Antibody, designed for immunoprecipitation (IP) of tagged protein
ABT2053	Agarose Conjugated Anti-His Tag Mouse Monoclonal Antibody (5C3)	Featured His Tag Mouse Monoclonal Antibody, designed for immunoprecipitation (IP) of tagged protein

## Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.