



FRS2 (phospho Tyr436) Polyclonal Antibody

Cat #: ABP50497

Size: 30µl /100µl /200µl

Product Information



	Product Name: FRS2 (phospho Tyr436) Polyclonal Antibody		
	Applications: WB , IHC-P, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse, Monkey		
REF	Catalog Number: ABP50497	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: Fibroblast growth factor (FGF) receptor substrate 2 (FRS2) has an alternative name as SNT-1, it is an adapter protein that links activated FGR and NGF receptors to downstream signaling pathways. FGF receptor substrates (FRS2 and FRS3) are key adaptor proteins that mediate FGF-FGFR signalling in benign as well as malignant tissue. FRS2 is a 508 amino-acid protein, which is phosphorylated on tyrosine residues. The molecular weight of non-phosphorylated FRS2 is 57-68 kDa, but phosphorylated FRS2 is 80-90 kDa. Phosphorylation of FRS2 is associated with activation of a number of MAP kinases. Allele-specific regulation of FGFR2 mRNA expression with a mildly increased breast cancer risk has been reported.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), IHC-P (1:100-1:300), ELISA (1:10000). Not yet tested in other applications.

Storage Buffer: PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

 <table data-bbox="325 351 651 441"><tr><td>-</td><td>+</td><td>- phospho-peptide</td></tr><tr><td>-</td><td>-</td><td>+ non-phospho-peptide</td></tr><tr><td>+</td><td>+</td><td>+ 3T3 NGF(customer's sample)</td></tr></table>	-	+	- phospho-peptide	-	-	+ non-phospho-peptide	+	+	+ 3T3 NGF(customer's sample)	<p>Fig.1. Western Blot analysis of various cells using Phospho-FRS2 (Y436) Polyclonal Antibody.</p>
-	+	- phospho-peptide								
-	-	+ non-phospho-peptide								
+	+	+ 3T3 NGF(customer's sample)								
	<p>Fig.2. Western Blot analysis of COS cells using Phospho-FRS2 (Y436) Polyclonal Antibody.</p>									

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.