

## TIRAP Polyclonal Antibody

Cat #: ABP53646

Size: 30µl /100µl /200µl

### Product Information

|                                                                                   |                                                                      |                                                                                   |                                           |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------|
|                                                                                   | <b>Product Name:</b> TIRAP Polyclonal Antibody                       |                                                                                   |                                           |
|                                                                                   | <b>Applications:</b> IHC-P, ELISA                                    |                                                                                   | <b>Isotype:</b> Rabbit IgG                |
|                                                                                   | <b>Reactivity:</b> Human, Mouse                                      |                                                                                   |                                           |
| <b>REF</b>                                                                        | <b>Catalog Number:</b> ABP53646                                      | <b>LOT</b>                                                                        | <b>Lot Number:</b> Refer to product label |
|                                                                                   | <b>Formulation:</b> Liquid                                           |                                                                                   | <b>Concentration:</b> 1 mg/ml             |
|  | <b>Storage:</b> Store at -20°C. Avoid repeated freeze / thaw cycles. |  | <b>Note:</b> Contain sodium azide.        |

**Background:** The innate immune system recognizes microbial pathogens through Toll-like receptors (TLRs), which identify pathogen-associated molecular patterns. Different TLRs recognize different pathogen-associated molecular patterns and all TLRs have a Toll-interleukin 1 receptor (TIR) domain, which is responsible for signal transduction. The TIR domain containing adaptor protein encoded by TIRAP is a TIR adaptor protein involved in the TLR4 signaling pathway of the immune system. It activates NF-kappa-B, MAPK1, MAPK3 and JNK, which then results in cytokine secretion and the inflammatory response. Alternative splicing of this gene results in several transcript variants; however, not all variants have been fully described.

**Application Notes:** Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IHC-P (1:100-1:300), ELISA (1:20000). Not yet tested in other applications.

**Storage Buffer:** PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.

**Storage Instructions:** Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

**Note:** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.