

## CaMKIV Polyclonal Antibody

Cat #: ABP56637

Size: 30µl /100µl /200µl

### Product Information

|                                                                                   |                                                                      |                                                                                   |                                           |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------|
|                                                                                   | <b>Product Name:</b> CaMKIV Polyclonal Antibody                      |                                                                                   |                                           |
|                                                                                   | <b>Applications:</b> WB, IHC-P, IF, ELISA                            |                                                                                   | <b>Isotype:</b> Rabbit IgG                |
|                                                                                   | <b>Reactivity:</b> Human, Mouse, Rat                                 |                                                                                   |                                           |
| <b>REF</b>                                                                        | <b>Catalog Number:</b> ABP56637                                      | <b>LOT</b>                                                                        | <b>Lot Number:</b> Refer to product label |
|                                                                                   | <b>Formulation:</b> Liquid                                           |                                                                                   | <b>Concentration:</b> 1 mg/ml             |
|  | <b>Storage:</b> Store at -20°C. Avoid repeated freeze / thaw cycles. |  | <b>Note:</b> Contain sodium azide.        |

**Background:** The product of CAMK4 (calcium/calmodulin dependent protein kinase IV) belongs to the serine/threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. This enzyme is a multifunctional serine/threonine protein kinase with limited tissue distribution, that has been implicated in transcriptional regulation in lymphocytes, neurons and male germ cells.

**Application Notes:** Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), IHC-P (1:100-1:300), IF (1:200-1:1000), ELISA (1:5000). Not yet tested in other applications.

**Storage Buffer:** PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.

**Storage Instructions:** Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

**Note:** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.