

CHD8 Polyclonal Antibody

Cat #: ABP58139

Size: 30µl /100µl /200µl

Product Information

	Product Name: CHD8 Polyclonal Antibody		
	Applications: IHC-P		Isotype: Rabbit IgG
	Reactivity: Human		
REF	Catalog Number: ABP58139	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: CHD8 (Chromodomain Helicase DNA Binding Protein 8) is a Protein Coding gene. Diseases associated with CHD8 include Autism 18 and Autism Spectrum Disorder Susceptibility, Chd8-Related. Among its related pathways are Regulation of Wnt-mediated beta catenin signaling and target gene transcription and Signaling by Wnt.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IHC-P (1:50-1:300).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.