

ENOF1 Polyclonal Antibody

Cat #: ABP58484

Size: 30µl /100µl /200µl

Product Information

	Product Name: ENOF1 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human		
REF	Catalog Number: ABP58484	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: ENOSF1 can encode a mitochondrial enzyme that is thought to convert L-fuconate to 2-keto-3-deoxy-L-fuconate. This locus was originally identified as the source of antisense RNAs of the adjacent thymidylate synthase gene. Splice variants at this locus may contain an alternate 3' exon that is complementary to the 3'UTR and terminal intron of the thymidylate synthase (TS) RNA and may downregulate TS expression. ENOSF1 (Enolase Superfamily Member 1) is a Protein Coding gene. Among its related pathways are Fructose and mannose metabolism.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.