

GDS1 Polyclonal Antibody

Cat #: ABP58628

Size: 30µl /100µl /200µl

Product Information

	Product Name: GDS1 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human		
REF	Catalog Number: ABP58628	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: RAP1GDS1 (Rap1 GTPase-GDP Dissociation Stimulator 1) is a Protein Coding gene. Diseases associated with RAP1GDS1 include Acute Lymphocytic Leukemia. Among its related pathways are G-protein signaling_Regulation of RAC1 activity and G-protein signaling_Rap1A regulation pathway. GO annotations related to RAP1GDS1 include binding and proton-transporting ATPase activity, rotational mechanism. The smg GDP dissociation stimulator (smgGDS) protein is a stimulatory GDP/GTP exchange protein with GTPase activity.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.