

K1C12 Polyclonal Antibody

Cat #: ABP58992

Size: 30µl /100µl /200µl

Product Information

	Product Name: K1C12 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse, Rat		
REF	Catalog Number: ABP58992	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: KRT12 (Keratin 12) is a Protein Coding gene. Diseases associated with KRT12 include Meesmann Corneal Dystrophy and Corneal Dystrophy. Among its related pathways are Keratinization and Developmental Biology. KRT12 encodes the type I intermediate filament chain keratin 12, expressed in corneal epithelia. Mutations in KRT12 lead to Meesmann corneal dystrophy.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.