



MYO15 Polyclonal Antibody

Cat #: ABP59373

Size: 30µl /100µl /200µl

Product Information

	Product Name: MYO15 Polyclonal Antibody		
	Applications: IHC-P		Isotype: Rabbit IgG
	Reactivity: Human, Mouse		
REF	Catalog Number: ABP59373	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: MYO15A (Myosin XVA) is a Protein Coding gene. Diseases associated with MYO15A include Deafness, Autosomal Recessive 3 and Nonsyndromic Deafness. Among its related pathways are Sertoli-Sertoli Cell Junction Dynamics and Actin Nucleation by ARP-WASP Complex. MYO15A encodes an unconventional myosin. This protein differs from other myosins in that it has a long N-terminal extension preceding the conserved motor domain.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IHC-P (1:50-1:300).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.