

## PO3F3 Polyclonal Antibody

Cat #: ABP59963

Size: 30µl /100µl /200µl

### Product Information

|                                                                                   |                                                                      |                                                                                   |                                           |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------|
|                                                                                   | <b>Product Name:</b> PO3F3 Polyclonal Antibody                       |                                                                                   |                                           |
|                                                                                   | <b>Applications:</b> WB, ELISA                                       |                                                                                   | <b>Isotype:</b> Rabbit IgG                |
|                                                                                   | <b>Reactivity:</b> Human, Mouse, Rat                                 |                                                                                   |                                           |
| <b>REF</b>                                                                        | <b>Catalog Number:</b> ABP59963                                      | <b>LOT</b>                                                                        | <b>Lot Number:</b> Refer to product label |
|                                                                                   | <b>Formulation:</b> Liquid                                           |                                                                                   | <b>Concentration:</b> 1 mg/ml             |
|  | <b>Storage:</b> Store at -20°C. Avoid repeated freeze / thaw cycles. |  | <b>Note:</b> Contain sodium azide.        |

**Background:** POU3F3 (POU Class 3 Homeobox 3) is a Protein Coding gene. Diseases associated with POU3F3 include Central Nervous System Tuberculosis and Agraphia. POU3F3 encodes a POU-domain containing protein that functions as a transcription factor. The encoded protein recognizes an octamer sequence in the DNA of target genes. This protein may play a role in development of the nervous system.

**Application Notes:** Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

**Storage Buffer:** PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

**Storage Instructions:** Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

**Note:** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.