

## RAC2 Polyclonal Antibody

Cat #: ABP60076

Size: 30µl /100µl /200µl

### Product Information

|                                                                                   |                                                                      |                                                                                   |                                           |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------|
|                                                                                   | <b>Product Name:</b> RAC2 Polyclonal Antibody                        |                                                                                   |                                           |
|                                                                                   | <b>Applications:</b> WB, ELISA                                       |                                                                                   | <b>Isotype:</b> Rabbit IgG                |
|                                                                                   | <b>Reactivity:</b> Human, Mouse                                      |                                                                                   |                                           |
| <b>REF</b>                                                                        | <b>Catalog Number:</b> ABP60076                                      | <b>LOT</b>                                                                        | <b>Lot Number:</b> Refer to product label |
|                                                                                   | <b>Formulation:</b> Liquid                                           |                                                                                   | <b>Concentration:</b> 1 mg/ml             |
|  | <b>Storage:</b> Store at -20°C. Avoid repeated freeze / thaw cycles. |  | <b>Note:</b> Contain sodium azide.        |

**Background:** RAC2 (Rac Family Small GTPase 2) encodes a member of the Ras superfamily of small guanosine triphosphate (GTP)-metabolizing proteins. The encoded protein localizes to the plasma membrane, where it regulates diverse processes, such as secretion, phagocytosis, and cell polarization. Activity of this protein is also involved in the generation of reactive oxygen species. Mutations in RAC2 are associated with neutrophil immunodeficiency syndrome. There is a pseudogene for RAC2 on chromosome 6.

**Application Notes:** Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

**Storage Buffer:** PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

**Storage Instructions:** Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

**Note:** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.