



RGL2 Polyclonal Antibody

Cat #: ABP60146

Size: 30µl /100µl /200µl

Product Information

	Product Name: RGL2 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse		
REF	Catalog Number: ABP60146	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: RGL2 (Ral Guanine Nucleotide Dissociation Stimulator Like 2) is a Protein Coding gene. Diseases associated with RGL2 include Tumoral Calcinosis, Familial, Normophosphatemic and Familial Tumoral Calcinosis. Among its related pathways are Ras signaling pathway and G-protein signaling N-RAS regulation pathway. GO annotations related to RGL2 include guanyl-nucleotide exchange factor activity and Ras guanyl-nucleotide exchange factor activity. An important paralog of RGL2 is RALGDS.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.