



RINT1 Polyclonal Antibody

Cat #: ABP60177

Size: 30µl /100µl /200µl

Product Information

	Product Name: RINT1 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse		
REF	Catalog Number: ABP60177	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: RINT1 (RAD50 Interactor 1) encodes a protein first identified for its ability to interact with the RAD50 double strand break repair protein, with the resulting interaction implicated in the regulation of cell cycle progression and telomere length. The encoded protein may also play a role in trafficking of cellular cargo from the endosome to the trans-Golgi network. Mutations in RINT1 may be associated with breast cancer in human patients.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.