



RIT2 Polyclonal Antibody

Cat #: ABP60182

Size: 30µl /100µl /200µl

Product Information

	Product Name: RIT2 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse, Rat		
REF	Catalog Number: ABP60182	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: RIN belongs to the RAS (HRAS; MIM 190020) superfamily of small GTPases (Shao et al. 1999 [PubMed 10545207]). RIT2 (Ras Like Without CAAX 2) is a Protein Coding gene. Diseases associated with RIT2 include Ehrlichiosis and Parkinson Disease Susceptibility. Among its related pathways are RET signaling and Neurotrophic factor-mediated Trk receptor signaling.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.