



RPA1 Polyclonal Antibody

Cat #: ABP60232

Size: 30µl /100µl /200µl

Product Information

	Product Name: RPA1 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human		
REF	Catalog Number: ABP60232	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: POLR1A (RNA Polymerase I Subunit A) is a Protein Coding gene. Diseases associated with POLR1A include Acrofacial Dysostosis, Cincinnati Type and Acrofacial Dysostosis. Among its related pathways are Mismatch repair and Pyrimidine metabolism (KEGG). The protein encoded by this gene is the largest subunit of the RNA polymerase I complex. The encoded protein represents the catalytic subunit of the complex, which transcribes DNA into ribosomal RNA precursors. Defects in POLR1A are a cause of the Cincinnati type of acrofacial dysostosis.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.