



## SFPQ Polyclonal Antibody

Cat #: ABP60374

Size: 30µl /100µl /200µl

### Product Information

	<b>Product Name:</b> SFPQ Polyclonal Antibody		
	<b>Applications:</b> WB, ELISA		<b>Isotype:</b> Rabbit IgG
	<b>Reactivity:</b> Human		
<b>REF</b>	<b>Catalog Number:</b> ABP60374	<b>LOT</b>	<b>Lot Number:</b> Refer to product label
	<b>Formulation:</b> Liquid		<b>Concentration:</b> 1 mg/ml
	<b>Storage:</b> Store at -20°C. Avoid repeated freeze / thaw cycles.		<b>Note:</b> Contain sodium azide.

**Background:** SFPQ (Splicing Factor Proline And Glutamine Rich) is a Protein Coding gene. Diseases associated with SFPQ include Renal Cell Carcinoma, Papillary and Adrenal Neuroblastoma. Among its related pathways are mRNA Splicing - Major Pathway and Circadian rhythm related genes. GO annotations related to this gene include nucleic acid binding and nucleotide binding. An important paralog of this gene is PSPC1.

**Application Notes:** Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

**Storage Buffer:** PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

**Storage Instructions:** Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

**Note:** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.