



## TAAR2 Polyclonal Antibody

Cat #: ABP60600

Size: 30µl /100µl /200µl

### Product Information

	<b>Product Name:</b> TAAR2 Polyclonal Antibody		
	<b>Applications:</b> WB, ELISA		<b>Isotype:</b> Rabbit IgG
	<b>Reactivity:</b> Human, Mouse, Rat		
<b>REF</b>	<b>Catalog Number:</b> ABP60600	<b>LOT</b>	<b>Lot Number:</b> Refer to product label
	<b>Formulation:</b> Liquid		<b>Concentration:</b> 1 mg/ml
	<b>Storage:</b> Store at -20°C. Avoid repeated freeze / thaw cycles.		<b>Note:</b> Contain sodium azide.

**Background:** TAAR2 (Trace Amine Associated Receptor 2 (Gene/Pseudogene)) is a Protein Coding gene. Diseases associated with TAAR2 include Generalized Epilepsy With Febrile Seizures Plus. Among its related pathways are G alpha(s) signalling events and Signaling by GPCR. GO annotations related to this gene include G-protein coupled receptor activity and trace-amine receptor activity. An important paralog of this gene is TAAR1.

**Application Notes:** Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

**Storage Buffer:** PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

**Storage Instructions:** Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

**Note:** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.