

TLX3 Polyclonal Antibody

Cat #: ABP60702

Size: 30µl /100µl /200µl

Product Information

	Product Name: TLX3 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse		
REF	Catalog Number: ABP60702	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: TLX3 (T-Cell Leukemia Homeobox 3) is a Protein Coding gene. Diseases associated with TLX3 include Precursor T-Cell Acute Lymphoblastic Leukemia and T-Cell Leukemia. Among its related pathways are Transcriptional misregulation in cancer and Sudden Infant Death Syndrome (SIDS) Susceptibility Pathways. The protein encoded by TLX3 is an orphan homeobox protein that encodes a DNA-binding nuclear transcription factor. A translocation [t(5;14)(q35;q32)] involving TLX3 is associated with T-cell acute lymphoblastic leukemia (T-ALL) in children and young adults.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.