

TPPC4 Polyclonal Antibody

Cat #: ABP60741

Size: 30µl /100µl /200µl

Product Information

	Product Name: TPPC4 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse, Rat		
REF	Catalog Number: ABP60741	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: TRAPPC4 (Trafficking Protein Particle Complex 4) is a Protein Coding gene. Among its related pathways are Syndecan-2-mediated signaling events and Transport to the Golgi and subsequent modification. GO annotations related to this gene include Rab guanyl-nucleotide exchange factor activity. May play a role in vesicular transport from endoplasmic reticulum to Golgi.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.