



TXLNA Polyclonal Antibody

Cat #: ABP60797

Size: 30µl /100µl /200µl

Product Information

	Product Name: TXLNA Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse		
REF	Catalog Number: ABP60797	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: TXLNA (Taxilin Alpha) is a Protein Coding gene. Diseases associated with TXLNA include Inclusion Conjunctivitis and Amyotrophic Lateral Sclerosis Type 10. Among its related pathways are Cytokine Signaling in Immune system and Innate Immune System. GO annotations related to this gene include cytokine activity and high molecular weight B cell growth factor receptor binding. An important paralog of this gene is TXLNG.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.