

UB2E2 Polyclonal Antibody

Cat #: ABP60809

Size: 30µl /100µl /200µl

Product Information

	Product Name: UB2E2 Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse		
REF	Catalog Number: ABP60809	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: UBE2E2 (Ubiquitin Conjugating Enzyme E2 E2) is a Protein Coding gene. Among its related pathways are Class I MHC mediated antigen processing and presentation and Ubiquitin mediated proteolysis. GO annotations related to this gene include ligase activity and acid-amino acid ligase activity. Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes Lys-11- and Lys-48-, as well as Lys-63-linked polyubiquitination. Catalyzes the ISGylation of influenza A virus NS1 protein.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.