

Fluo-3, AM ester

Cat #:BMD00006

Lot #: Refer to product label

Membrane-permeant AM ester form of Fluo-3 that can be loaded into cells via incubation.

Product Information

Product name Fluo-3, AM ester

Specification




Fig. Fluo-3, AM ester structure formula

Product Properties

Storage instructions

Store at -20°C. The product can be stored for up to 12 months. Protect from light.

Shipping

Gel pack with blue ice.

Precautions

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Additional Information

Background

Membrane-permeant AM ester form of Fluo-3 that can be loaded into cells via incubation. Fluo-3 AM ester itself does not bind Ca²⁺, but it is readily hydrolyzed to fluo-3 by endogenous esterases once the dye is inside the cells.

For more latest publications, please visit website www.abbkine.com