

6-TAMRA

Cat #:BMD00032

Lot #: Refer to product label

Tetramethylrhodamine has been a widely used fluorophore for preparing bioconjugates.

Product Information

Product name 6-TAMRA

Specification




Fig. 6-TAMRA structure
formula

Product Properties

Storage instructions	Store at +4°C. The product can be stored for up to 12 months. Protect from light.
Shipping	Gel pack with blue ice.
Precautions	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Additional Information

Background	Tetramethylrhodamine (TMR, TRITC) has been a widely used fluorophore for preparing bioconjugates, such as peptide, protein, nucleotide and nucleic acid conjugation, especially fluorescent antibody and avidin derivatives used in immunochemistry. 6-TAMRA is one of the traditional fluorophores used for automated DNA sequencing.
-------------------	--

For more latest publications, please visit website www.abbkine.com