

## 5-ROX

Cat #:BMD00035

Lot #: Refer to product label

5-ROX is used for oligonucleotide labeling and automated DNA sequencing applications.

### Product Information

**Product name** 5-ROX

### Specification




Fig. 5-ROX structure formula

### Product Properties

**Storage instructions** Store at +4°C. The product can be stored for up to 12 months. Protect from light.

**Shipping** Gel pack with blue ice.

**Precautions** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

### Additional Information

**Background** The carboxylic acid of 5-ROX is used for oligonucleotide labeling and automated DNA sequencing applications.

**Alternative names** 5-carboxy-X-rhodamine

For more latest publications, please visit website [www.abbkine.com](http://www.abbkine.com)