

6-ROX

Cat #:BMD00036

Lot #: Refer to product label

6-ROX is predominately used for labeling nucleotides and sequencing nucleic acids.

Product Information

Product name 6-ROX

Specification




Fig. 6-ROX structure formula

Product Properties

Storage instructions	Store at -20°C. The product can be stored for up to 12 months. Protect from light.
Shipping	Gel pack with blue ice.
Precautions	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Additional Information

Background	6-ROX is predominately used for labeling nucleotides and sequencing nucleic acids.
Alternative names	6-Carboxy-X-rhodamine

For more latest publications, please visit website www.abbkine.com