

6-ROX, succinimidyl ester

Cat #:BMD00037

Lot #: Refer to product label

6-ROX, SE is widely used for oligonucleotide labeling and automated DNA sequencing applications.

Product Information

Product name 6-ROX, succinimidyl ester

Specification




Fig. 6-ROX, succinimidyl ester structure formula

Product Properties

Storage instructions	Store at -20°C. The product can be stored for up to 12 months. Protect from light.
Shipping	Gel pack with blue ice.
Precautions	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Additional Information

Background	The isomeric 6-ROX, SE is an amine-reactive form of carboxy-X-rhodamine and is widely used for oligonucleotide labeling and automated DNA sequencing applications.
-------------------	--

For more latest publications, please visit website www.abbkine.com