

## 5-FAM, succinimidyl ester

Cat #:BMD00043

Lot #: Refer to product label

5-FAM, SE is the amine-reactive succinimidyl ester of single isomer 5-FAM acid.

### Product Information

**Product name** 5-FAM, succinimidyl ester

### Specification




Fig. 5-FAM, succinimidyl ester structure formula

### Product Properties

**Storage instructions** Store at -20°C. The product can be stored for up to 12 months. Protect from light.

**Shipping** Gel pack with blue ice.

**Precautions** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

### Additional Information

**Background** 5-FAM, SE is the amine-reactive succinimidyl ester of single isomer 5-FAM acid. It is one of the most popular green fluorescent reagents used for labeling peptides, proteins and nucleotides. It has also been used to prepare various small fluorescent molecules.

**Alternative names** 5-Carboxyfluorescein, Succinimidyl Ester

For more latest publications, please visit website [www.abbkine.com](http://www.abbkine.com)