

5-FAM-PEO8, succinimidyl ester

Cat #:BMD00044

Lot #: Refer to product label

5-FAM-PEO8 SE has an approximately 32 Angstrom water-soluble PEG spacer that separates fluorescein and the amine-reactive SE group.

Product Information

Product name 5-FAM-PEO8, succinimidyl ester

Specification




Fig. 5-FAM-PEO8, succinimidyl ester structure formula

Product Properties

Storage instructions	Store at -20°C. The product can be stored for up to 12 months. Protect from light.
Shipping	Gel pack with blue ice.
Precautions	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Additional Information

Background	5-FAM-PEO8 SE (full name: 5-Carboxyfluorescein-PEO8 propionate succinimidyl ester) has an approximately 32 Angstrom water-soluble PEG spacer that separates fluorescein and the amine-reactive SE group.
-------------------	--

For more latest publications, please visit website www.abbkine.com