

6-FAM, succinimidyl ester

Cat #:BMD00047

Lot #: Refer to product label

6-FAM, SE is the amine-reactive succinimidyl ester of single isomer 6-FAM acid.

Product Information

Product name 6-FAM, succinimidyl ester

Specification




Fig. 6-FAM, succinimidyl ester structure formula

Product Properties

Storage instructions	Store at -20°C. The product can be stored for up to 12 months. Protect from light.
Shipping	Gel pack with blue ice.
Precautions	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Additional Information

Background	6-FAM, SE is the amine-reactive succinimidyl ester of single isomer 6-FAM acid. It is one of the most popular green fluorescent reagents used for labeling nucleotides and nucleic acids. Compared to 5-FAM, 6-FAM is less often used to prepare small molecules.
Alternative names	6-Carboxyfluorescein, Succinimidyl Ester

For more latest publications, please visit website www.abbkine.com