

Human CNTF protein

Cat #: PRP1002

Size: 5 µg/20 µg/100 µg/1 mg

Product Information

	Product Name: Human CNTF protein		
REF	Catalog Number: PRP1002	LOT	Lot Number: Refer to product label
	Purity: >98% as determined by SDS-PAGE		
	Storage: Store at -20°C		Preparation method: <i>E.Coli</i>
	Shipping: The product is shipped at ambient temperature.		

Background: Ciliary neurotrophic factor (CNTF) is a member of the cytokine family. CNTF is a survival factor for various neuronal cell types. Subsequent studies have demonstrated that CNTF is a survival factor for additional neuronal cell types including: dorsal root ganglion sensory neurons, sympathetic ganglion neurons, embryonic motor neurons, major pelvic ganglion neurons and hippocampal neurons.

Sequence: Amino acid sequence derived from Human CNTF(P26441)(Ala2-Met200) was expressed.

Protein length: The recombinant Human CNTF consists of 199 amino acids and has a predicted molecular mass of 22 kDa.


Biological Activity: Test in progress.

Endotoxin: <0.1 EU per µg of the protein as determined by the LAL method.

Formulation: Lyophilized from sterile PBS, pH 7.5.

Storage Instructions: Lyophilized Human CNTF protein product should be stored desiccated below -18°C. Upon reconstitution, the protein should be stored at 4°C between 2 -7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage notes: Always centrifuge tubes before opening. It is recommended to reconstitute the lyophilized Human CNTF protein using the buffer we provided not less than 100 µg/mL, which can then be further diluted to other aqueous solutions.



Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We can not be responsible for patent infringements or other violations that may occur with the use of this product.