


Mouse CNTF Protein, His tag (Animal-Free)

Cat #: PRP1140

Size: 5µg/20µg/100µg/1mg

Product Information

	Product Name: Mouse CNTF Protein, His tag (Animal-Free)		
REF	Catalog Number: PRP1140	LOT	Lot Number: Refer to product label
	Purity: > 98% as determined by SDS-PAGE		
	Storage: Store at -20°C		Preparation method: E.coli
	Shipping: The product is shipped at ambient temperature.		

Background: Ciliary neurotrophic factor (CNTF) is a member of the cytokine family. It is a polypeptide hormone that has functions in promoting neurotransmitter synthesis and neurite outgrowth in certain neuronal populations. Its actions appear to be restricted to the nervous system. Ciliary neurotrophic factor (CNTF) has biological effects through the activation of a multi-subunit receptor complex, consisting of an extracellular CNTF binding subunit (CNTF α) and two transmembrane signal transduction proteins: glycoprotein gp130 and LIF receptor. CNTF is considered as a potent survival factor of neurons and oligodendrocyte and may be relevant in reducing tissue destruction during inflammatory attacks.

Sequence: Amino acid sequence derived from Mouse CNTF (Ala2-Met198)(P51642) was expressed with 6×His tag at the C-terminus.

Protein length: The recombinant Mouse CNTF consists of 198 amino acids and predicts a molecular mass of 22.6 kDa.

Biological Activity: Testing in progress.


Endotoxin: < 0.1 EU per µg of the protein as determined by the LAL method.

Formulation: The protein was lyophilized from a 0.2 µm filtered solution containing 1X PBS, pH 7.4. If you have any concerns or special requirements, please confirm with us.

Storage Instructions: Lyophilized protein product should be stored desiccated below -18°C. Upon reconstitution, the protein should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (5% HSA, 10% FBS or 0.1% BSA). Please prevent freeze-thaw cycles.

Usage notes: Always centrifuge tubes before opening. It is recommended to reconstitute the lyophilized recombinant Mouse

CNTF protein to a concentration of 0.1-1 mg/mL in sterile ddH₂O, and keep at room temperature for at least 20 min to fully dissolve. Please avoid vortex vigorously.



Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.