


Rat IL-6 Protein

Cat #: PRP1241

Size: 5µg/20µg/100µg/1mg

Product Information

	Product Name: Rat IL-6 Protein		
REF	Catalog Number: PRP1241	LOT	Lot Number: Refer to product label
	Purity: > 99 % as determined by SDS-PAGE		
	Storage: Store at -20°C		Preparation method: <i>E. coli</i>
	Shipping: The product is shipped at ambient temperature.		

Background: Interleukin-6 (IL-6) is a multifunctional α -helical cytokine that regulates cell growth and differentiation of various tissues, which is known particularly for its role in the immune response and acute phase reactions. IL-6 protein is secreted by a variety of cell types including T cells and macrophages as a phosphorylated and variably glycosylated molecule. It exerts actions through its heterodimeric receptor composed of IL-6R that lacks the tyrosine/kinase domain and binds IL-6 with low affinity, and ubiquitously expressed glycoprotein 130 (gp130) that binds the IL-6. IL-6R complex with high affinity and thus transduces signals. IL-6 is also involved in hematopoiesis, bone metabolism, and cancer progression, and has been defined as an essential role in directing the transition from innate to acquired immunity.

Sequence: Amino acid sequence derived from Recombinant Rat IL-6 Protein P20607 (F25-T211) was expressed.

Protein length: The Recombinant Rat IL-6 Protein consists of 187 amino acids and predicts a molecular mass of 21.73KDa


Biological Activity: Testing in progress.

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Formulation: Lyophilized from sterile 20 mM Tris with 150 mM NaCl, pH 7.4.

Storage Instructions: Lyophilized Recombinant Rat IL-6 Protein product should be stored desiccated below -20°C. Upon reconstitution, the protein should be stored at 4°C between 2-7 days and for future use below -20°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage notes: Always centrifuge tubes before opening. It is recommended to reconstitute the lyophilized Recombinant Rat IL-6 Protein using the buffer we provided not less than 100µg/ml, which can then be further diluted to other aqueous solutions.



Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.