


## Human CXCL9 Protein, His tag (Animal-Free)

Cat #: PRP1915

Size: 5 µg/20 µg/100 µg/1 mg

### Product Information

	<b>Product Name:</b> Human CXCL9 Protein, His tag (Animal-Free)		
<b>REF</b>	<b>Catalog Number:</b> PRP1915	<b>LOT</b>	<b>Lot Number:</b> Refer to product label
	<b>Purity:</b> > 98 % as determined by SDS-PAGE		
	<b>Storage:</b> Store at -20°C		<b>Preparation method:</b> E.coli
	<b>Shipping:</b> The product is shipped at ambient temperature.		

**Background:** C-X-C motif chemokine 9 (CXCL9) also named monokine induced by gamma interferon (MIG), which is a chemokine of the intercrine alpha family. CXCL9 controls the immune cells by binding the CXCR3 which is including the cell migration and activation. During inflammation, CXCL9 is a chemotaxis for lymphocyte and macrophages. CXCL9 is participated in the process of tumor proliferation and metastasis.

**Sequence:** Amino acid sequence derived from Human CXCL9 (Thr23-Thr125) (Q07325) was expressed with 6×His tag at the N-terminus.

**Protein length:** The recombinant Human CXCL9 consists of 103 amino acids and predicts a molecular mass of 12.5 kDa.


**Biological Activity:** Testing in process.

**Endotoxin:** <0.1 EU per 1 µg of the protein by the LAL method.

**Formulation:** Lyophilized from PBS, pH 7.4.

**Storage Instructions:** Lyophilized protein product should be stored desiccated below -18°C. Upon reconstitution, the protein should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (5% HSA, 10%FBS or 0.1%BSA). Please prevent freeze-thaw cycles.

**Usage notes:** Always centrifuge tubes before opening. It is recommended to reconstitute the lyophilized recombinant protein to a concentration of 0.1-1 mg/mL using the buffer we provided, and keep at room temperature for at least 20 min to fully dissolve. Please avoid vortex vigorously.



**Note:** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.