


Human CD73/NT5E protein, His tag (Animal-Free)

Cat #: PRP2036

Size: 5 µg/ 20 µg/100 µg/1 mg

Product Information

	Product Name: Human CD73 protein (Animal-Free)		
REF	Catalog Number: PRP2036	LOT	Lot Number: Refer to product label
	Purity: >95% as determined by SDS-PAGE		
	Storage: Store at -20C		Preparation method: HEK 293 Cells
	Shipping: The product is shipped at ambient temperature.		

Background: CD73 (ecto-5'-nucleotidase) is a novel immunoinhibitory protein that plays a key role for tumor growth and metastasis. Its main function is to convert extracellular ATP to immunosuppressive adenosine in concert with CD39 in normal tissues to limit excessive immune response. However, tumors take advantage of the CD73-mediated adenosinergic mechanism to protect them from immune attack. In particular, inducible expression of CD73 along with other adenosinergic molecules on both cancer cells and host cells sustains immunosuppressive tumor microenvironment by affecting multiple aspects of the immune response. Owing to its multifaceted capacity to tumor promotion as an emerging immune checkpoint, CD73 is an ideal therapeutic target for cancer treatment especially in combination with conventional therapy and/or other immune checkpoint inhibitors.

Sequence: Amino acid sequence derived from human CD73 protein (NP_002517.1, Met1-Lys547) was expressed with a polyhistidine tag at the C-terminus.

Protein length: The recombinant human CD73 protein consists of 528 amino acids and predicts a molecular mass of 59 kD. The apparent molecular mass of the human CD73 is approximately 70 kD in SDS-PAGE under reducing conditions due to glycosylation.


Biological Activity: Testing in progress.

Endotoxin: < 1 EU per µg of the protein as determined by the LAL method.

Formulation: Lyophilized from sterile PBS, pH 7.4.

Storage Instructions: Lyophilized human CD73 protein product should be stored desiccated below -18°C. Upon reconstitution, the protein should be stored at 4°C between 2 -7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage notes: Always centrifuge tubes before opening. It is recommended to reconstitute the lyophilized human CD73 protein in the buffer which we provide not less than 100 µg/mL, which can then be further diluted to other aqueous solutions.



Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.