

Abbkine Top Selling Products for Cytology

Abbkine focuses on Proteomics and Cytology. Here, we present our top selling products in cytology research, from dyes and kits for cell status detection, organelle extraction kits, cell substructure staining and cell metabolism detection products, to cytokine for cell culture and protein detection kits, to help your research career. The following products are available in stock.

Cell Status

Cell Status Detection Kit (Proliferation & Toxicity & Apoptosis & Senescence)

Product Name	Description	Cat. No.	Size
Cell Counting Kit-8 (CCK-8)	Cell Proliferation/Toxicity	KTC011001	100/500/2000/10000T
LDH Cytotoxicity Assay Kit	Cell Toxicity	KTA1030	100/500/2000T
Cell Senescence Kit (β -Galactosidase Assay)	Cell Senescence	KTA3010	100/500T
TUNEL Apoptosis Detection Kit	Cell Apoptosis (DNA Fragmentation)	KTA2010	20/50/100T
Mitochondrial Membrane Potential Assay Kit (JC-1)	Cell Apoptosis (Mitochondrial Membrane Potential)	KTA4001	20/100/500T
Annexin V-AbFluor™ 405 Apoptosis Detection Kit	Cell Apoptosis (Blue Fluorescence)	KTA0001	20/50/100T
Annexin V-AbFluor™ 488 Apoptosis Detection Kit	Cell Apoptosis (Green Fluorescence)	KTA0002	20/50/100T
Annexin V-AbFluor™ 555 Apoptosis Detection Kit	Cell Apoptosis (Orange Fluorescence)	KTA0003	20/50/100T
Annexin V-AbFluor™ 647 Apoptosis Detection Kit	Cell Apoptosis (Deep Red Fluorescence)	KTA0004	20/50/100T
Annexin V-AbFluor™ 488 Reagent	Cell Apoptosis Reagent (Green Fluorescence)	BMD10082	200/1000T
Annexin V-AbFluor™ 555 Reagent	Cell Apoptosis Reagent (Orange Fluorescence)	BMD10083	200/1000T

*** Combine Annexin V with PI to detect early apoptotic and necrotic cells. **Note:** KTA0003 does not contain PI. More [cell status detection products](#), please contact Abbkine or local distributors.

Organelle Extraction

ExKine™ Organelle Extraction Kit (Complete & Active & Convenient)

Product Name	Description	Cat. No.	Size
ExKine™ Nuclei Extraction Kit	Nuclei Extraction (Crude)	KTP4001	50/200T
ExKine™ Nuclei Extraction Kit (High Purity)	Nuclei Extraction (High Purity)	KTP4002	20/100T
ExKine™ Mitochondrion Extraction Kit (Cultured Cells)	Mitochondrion Extraction (Cells)	KTP4003	50/200T
ExKine™ Mitochondrion Extraction Kit (Tissue)	Mitochondrion Extraction (Tissue)	KTP4004	50/200T
ExKine™ Nuclear and Cytoplasmic Protein Extraction Kit	Nuclear & Cytoplasmic Protein Extraction	KTP3001	50/200T
ExKine™ Membrane and Cytoplasmic Protein Extraction kit	Membrane & Cytoplasmic Protein Extraction	KTP3005	50/200T

*** Provide organelle and protein extraction kits to meet different needs, including nuclear protein, cytoplasmic protein, membrane protein extraction, etc. More [organelle or protein extraction products](#), please contact Abbkine or local distributors.

TraKine™ Cell Staining Kit (Low Toxicity & Anti-quenching & Co-staining)

Product Name	Description	Cat. No.	Size
Live and Dead Cell Double Staining Kit	Live and Dead Cell Double Staining	KTA1001	100/500/2000T
Live Cell Tracking Kit (Green Fluorescence)	Live Cell Tracking	KTA1002	100/500/2000T
TraKine™ Cell Plasma Membrane Staining Kit (Green Fluorescence)	Cell Plasma Membrane Staining	KTC4001	100/500/2000T
TraKine™ Cell Plasma Membrane Staining Kit (Orange Fluorescence)	Cell Plasma Membrane Staining	KTC4002	100/500/2000T
TraKine™ Mitochondrion Staining Kit (Green Fluorescence)	Mitochondrion Staining	KTC4003	100/500/2000T
TraKine™ Mitochondrion Staining Kit (Orange Fluorescence)	Mitochondrion Staining	KTC4004	100/500/2000T
TraKine™ Mitochondrion and Nuclear Staining Kit	Mitochondrion and Nuclear Staining	KTC4005	100/500/2000T
TraKine™ F-actin Staining Kit (Green Fluorescence)	F-actin Staining	KTC4008	50/300/1000T
TraKine™ F-actin Staining Kit (Orange Fluorescence)	F-actin Staining	KTC4009	50/300/1000T

*** Provide various cell and organelle staining and tracing kits. More [Cell Staining Products](#), please contact Abbkine or local distributors.

Fluorescent Dyes (Strong Light Resistance & Low Toxicity)

Product Name	Description	Cat. No.	Size
AbFluor™ 488-Phalloidin	Phalloidin (Green Fluorescence)	BMD00082	500T
AbFluor™ 555-Phalloidin	Phalloidin (Orange Fluorescence)	BMD00083	500T
AbFluor™ 594-Phalloidin	Phalloidin (Red Fluorescence)	BMD00084	500T
Hoechst 33342	Nuclei (Live Cell)	BMD00062	10/100mg
DAPI	Nuclei	BMD00063	10/100mg
Calcein AM	Cell Viability Detection (Live Cell)	BMD00064	1/10mg
6-CDCFDA, SE	Cell Viability Detection	BMD00070	5/50mg
DiI (DiIc18(3))	Cell Plasma Membrane (Orange Fluorescence)	BMD00071	20/200mg
DiO(DiOc18(3))	Cell Plasma Membrane (Green Fluorescence)	BMD00072	20/200mg
DiD (DiIc18(5))	Cell Plasma Membrane (Red Fluorescence)	BMD00073	10/100mg
DiR (DiIc18(7))	Cell Plasma Membrane (Deep Red Fluorescence)	BMD00074	5/50mg
DiOC6(3)	Organelle Membrane	BMD00081	20/200mg
RH421	Neuron Cell	BMD00076	5/50mg
Fluo-3, AM ester	Calcium Ion Fluorescent Probe	BMD00006	1/10mg
BAPTA, AM	Calcium Chelating Agent	BMD00085	25/250mg
MQAE	Chloride Ion Fluorescent Probe	BMD00008	50/500mg
BCECF, AM ester	pH Fluorescent Probe	BMD00007	1/10mg

*** SE groups can react with amino groups to produce stable amide bonds, which are more stable and better water soluble.

*** Provide fluorescent dyes and probes for various cells and organelles. More [fluorescent products](#), please contact Abbkine or local distributors.

CheKine™ Cell Metabolism Assay Kit (High Sensitivity & Stability)

Product Name	Description	Cat. No.	Size
CheKine™ NADP/NADPH Assay Kit	NADP/NADPH Quantification/Ratio Determination	KTB1010	48/96T
CheKine™ NAD/NADH Assay Kit	NAD/NADH Quantification/Ratio Determination	KTB1020	48/96T

*** More [cell metabolism assay products](#), please contact Abbkine or local distributors.

Cytokines (Low Endotoxin & High Purity & Freeze-dried Powder)

Product Name	Description	Cat. No.	Size
Human VEGF121 protein	VEGF121	PRP100005	5ug/20ug/100ug/1mg
Human aFGF protein	aFGF	PRP100008	10ug/50ug/100ug/1mg
Human bFGF protein	bFGF	PRP100009	10ug/50ug/100ug/1mg
Human GM-CSF protein	GM-CSF	PRP100013	5ug/20ug/100ug/1mg
Human IL-1 alpha protein	IL-1 α	PRP100048	5/20/100ug
Human TNF- β protein	TNF- β	PRP100078	5ug/20ug/100ug/1mg
Human IL-6 protein	IL-6	PRP100120	5 μ g/20 μ g/100 μ g/1mg
Human IGF1 protein	IGF1	PRP100153	100ug/200ug/500ug/1mg
Human TNF- α protein	TNF- α	PRP100157	10ug/50ug/100ug/1mg
Human EGF protein	EGF	PRP100159	100ug/500ug/1mg
Human IFN- β protein	IFN- β	PRP100183	5ug/20ug/100ug/1mg
Human TGFB1 protein	TGF- β 1	PRP100190	5ug/20ug/100ug/1mg
Human IL-10 protein	IL-10	PRP100228	5ug/20ug/100ug/1mg
Human RANKL protein	RANKL	PRP100262	5/20/100ug
Human IFN- γ protein	IFN- γ	PRP100263	20ug/100ug/500ug/1mg
Human IL-4 protein	IL-4	PRP100273	5ug/20ug/100ug/1mg
Human IL-17 protein	IL-17	PRP100290	5ug/20ug/100ug/1mg
Human IL-5 protein	IL-5	PRP100333	5/20/100ug
Mouse IL-6 protein	IL-6	PRP100371	5/20/100ug
Mouse IL-21 protein	IL-21	PRP100372	5/20/100ug
Mouse TNF- α protein	TNF- α	PRP100409	5ug/20ug/100ug/1mg
Mouse IL-4 protein	IL-4	PRP100498	5/20/100ug
Human IL-12 (IL-12A & IL-12B) protein	IL-12	PRP100705	5/20/50ug
Human Leptin protein	Leptin	PRP100730	5ug/500ug/1mg/2mg

*** Provide high-quality and highly biologically active cytokines. More [cytokine products](#), please contact Abbkine or local distributors.

LinKine™ Labeling Kit (Reliable & Efficiency)

Product Name	Description	Cat. No.	Size
LinKine™ HRP Labeling Kit	HRP Labeling	KTL0100	3x20ug/100ug/3x100ug/1mg
LinKine™ FITC Labeling Kit	FITC Labeling	KTL0210	3x20ug/100ug/3x100ug/1mg
LinKine™ AbFluor™ 350 Labeling Kit	AbFluor™350 (Blue Fluorescence) Labeling	KTL0500	3x20ug/100ug/3x100ug/1mg
LinKine™ AbFluor™ 405 Labeling Kit	AbFluor™405 (Blue Fluorescence) Labeling	KTL0510	3x20ug/100ug/3x100ug/1mg
LinKine™ AbFluor™ 488 Labeling Kit	AbFluor™488 (Green Fluorescence) Labeling	KTL0520	3x20ug/100ug/3x100ug/1mg
LinKine™ AbFluor™ 555 Labeling Kit	AbFluor™555 (Orange Fluorescence) Labeling	KTL0530	3x20ug/100ug/3x100ug/1mg
LinKine™ AbFluor™ 594 Labeling Kit	AbFluor™594 (Red Fluorescence) Labeling	KTL0540	3x20ug/100ug/3x100ug/1mg
LinKine™ AbFluor™ 647 Labeling Kit	AbFluor™647 (Deep Red Fluorescence) Labeling	KTL0560	3x20ug/100ug/3x100ug/1mg
LinKine™ AbFluor™ 680 Labeling Kit	AbFluor™680 (Near-infrared) Labeling	KTL0580	3x20ug/100ug/3x100ug/1mg
LinKine™ AbFluor™ 770 Labeling Kit	AbFluor™770 (Far-infrared) Labeling	KTL0590	3x20ug/100ug/3x100ug/1mg

*** Provide a full range of protein, antibody labeling schemes. More [protein couple products](#), please contact Abbkine or local distributors.

EliKine™ ELISA Kit (High Sensitivity & High Specificity)

Product Name	Detection Index	Cat. No.	Size
EliKine™ Human β -hCG ELISA Kit	Human β -hCG	KET6901	48/96/96x5T
EliKine™ Human LH ELISA Kit	Human LH	KET6902	48/96/96x5T
EliKine™ Human FSH ELISA Kit	Human FSH	KET6903	48/96/96x5T
EliKine™ Human Prolactin ELISA Kit	Human Prolactin	KET6904	48/96/96x5T
EliKine™ Human TSH ELISA Kit	Human TSH	KET6905	48/96/96x5T
EliKine™ Testosterone ELISA Kit	Testosterone	KET0001	48/96/96x5T
EliKine™ Progesterone ELISA Kit	Progesterone	KET0002	48/96/96x5T
EliKine™ Estradiol ELISA Kit	Estradiol	KET0003	48/96/96x5T
EliKine™ Free Triiodothyronine (FT3) ELISA Kit	Free Triiodothyronine	KET0004	48/96/96x5T
EliKine™ Free Thyroxine (FT4) ELISA Kit	Free Thyroxine	KET0005	48/96/96x5T
EliKine™ Triiodothyronine (T3) ELISA Kit	Triiodothyronine	KET0006	48/96/96x5T
EliKine™ Thyroxine (T4) ELISA Kit	ThyroxineC	KET0007	48/96/96x5T
EliKine™ Human CCL2 ELISA Kit	Human CCL2	KET6001	48/96/96x5T
EliKine™ Human EGF ELISA Kit	Human EGF	KET6005	48/96/96x5T
EliKine™ Human IL-1 α ELISA Kit	Human IL-1 α	KET6012	48/96/96x5T
EliKine™ Human IL-1 β ELISA Kit	Human IL-1 β	KET6013	48/96/96x5T
EliKine™ Human IL-2 ELISA Kit	Human IL-2	KET6014	48/96/96x5T
EliKine™ Human IL-6 ELISA Kit	Human IL-6	KET6017	48/96/96x5T
EliKine™ Human IL-8 ELISA Kit	Human IL-8	KET6018	48/96/96x5T
EliKine™ Human IL-10 ELISA Kit	Human IL-10	KET6019	48/96/96x5T
EliKine™ Human IL-17 ELISA Kit	Human IL-17	KET6022	48/96/96x5T
EliKine™ Human IL-18 ELISA Kit	Human IL-18	KET6023	48/96/96x5T
EliKine™ Human MMP-9 ELISA Kit	Human MMP-9	KET6027	48/96/96x5T
EliKine™ Human TGF- β 1 ELISA Kit	Human TGF- β 1	KET6030	48/96/96x5T
EliKine™ Human TIMP-1 ELISA Kit	Human TIMP-1	KET6031	48/96/96x5T
EliKine™ Human TNF- α ELISA Kit	Human TNF- α	KET6032	48/96/96x5T
EliKine™ Human VEGF ELISA Kit	Human VEGF	KET6033	48/96/96x5T
EliKine™ Mouse IFN- γ ELISA Kit	Mouse IFN- γ	KET7003	48/96/96x5T
EliKine™ Mouse IL-1 β ELISA Kit	Mouse IL-1 β	KET7005	48/96/96x5T
EliKine™ Mouse IL-2 ELISA Kit	Mouse IL-2	KET7006	48/96/96x5T
EliKine™ Mouse IL-6 ELISA Kit	Mouse IL-6	KET7009	48/96/96x5T
EliKine™ Mouse IL-10 ELISA Kit	Mouse IL-10	KET7010	48/96/96x5T
EliKine™ Mouse IL-12 p70 ELISA Kit	Mouse IL-12 p70	KET7011	48/96/96x5T
EliKine™ Mouse IL-17 ELISA Kit	Mouse IL-17	KET7012	48/96/96x5T
EliKine™ Mouse TGF- β 1 ELISA Kit	Mouse TGF- β 1	KET7014	48/96/96x5T
EliKine™ Mouse TNF- α ELISA Kit	Mouse TNF- α	KET7015	48/96/96x5T
EliKine™ Mouse VEGF ELISA Kit	Mouse VEGF	KET7016	48/96/96x5T
EliKine™ Rat IL-1 β ELISA Kit	Rat IL-1 β	KET9001	48/96/96x5T
EliKine™ Rat IL-6 ELISA Kit	Rat IL-6	KET9004	48/96/96x5T
EliKine™ Rat TNF- α ELISA Kit	Rat TNF- α	KET9007	48/96/96x5T

*** EliKine™ series of kits selected popular cytokines and hormones as detection indicators. More [ELISA kits](#), please contact Abbkine or local distributors.

To get more information, Please contact Abbkine or visit <https://www.abbkine.com>.

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.

