

Abbkine Top Selling Products for Immunology

Abbkine focuses on Immunology and Cytology. Here, we present our top selling products in immunology research, from basic Immunology products such as Protein Extraction and Quantification kits to Tag Antibodies, Primary Antibodies and Secondary Antibodies, to help your research career. The following products are available in stock.

Western Blot

Loading Control Antibodies

Product Name	Description	Cat. No.	Size
Anti-beta Actin Mab (1C7)	Universal Loading Control	A01010	50ul/200ul/1ml/10ml
HRP Conjugated Anti-beta Actin Mab (1C7)	HRP Conjugated	A01015	50ul/200ul/1ml
Anti-β-Actin Mab (1C7), AbFluor™ 488	Green Dye Conjugated	A01010A488	100ul
Anti-β-Actin Pab (A283)	WB, IHC	ABP57456	30ul/100ul/200ul
Anti-GAPDH Mab (2B5)	Universal Loading Control	A01020	50ul/200ul/1ml/10ml
Anti-GAPDH Pab	WB, IHC	ABP50152	30ul/100ul/200ul
HRP Conjugated Anti-GAPDH Mab (2B5)	HRP Conjugated	A01025	50ul/200ul/1ml
Anti-GAPDH Mab (2B5), AbFluor™ 488	Green Dye Conjugated	A01020A488	100ul
Anti-β-Tubulin Mab (3G6)	Universal Loading Control	A01030	50ul/200ul/1ml/10ml
Anti-β-Tubulin Mab (3G6), HRP Conjugate	HRP Conjugated	A01030HRP	50ul/200ul/1ml
Anti-β-Tubulin Mab (3G6), AbFluor™ 555	Orange Dye Conjugated	A01030A555	100ul
Anti-α-Tubulin Mouse Mab (3G5)	Universal Loading Control	A01080	50ul/200ul/1ml/10ml
Anti-PCNA Mab (1D7)	Nuclear Protein Loading Control	A01040	50ul/200ul/1ml/10ml
Anti-Histone H3 Mab (2D10)	Nuclear Protein Loading Control	A01070	50ul/200ul/1ml/10ml
Anti-Histone H3 Mab (2D10), AbFluor™ 488	Green Dye Conjugated	A01070A488	100ul
Anti-Lamin B1 Mab (15T1)	Nuclear Protein Loading Control	A01090	50ul/200ul/1ml/10ml
Anti-Lamin B1 Mab (15T1), AbFluor™ 594	Red Dye Conjugated	A01090A594	100ul
Anti-TBP/TATA Binding Protein Mab (2C6)	Nuclear Protein Loading Control	A01120	50ul/200ul/1ml/10ml
Anti-COX IV Mab (14Y2)	Mitochondrial Protein Loading Control	A01060	50ul/200ul/1ml/10ml
Anti-COX IV Mab (14Y2), HRP Conjugate	HRP Conjugated	A01060HRP	50ul/200ul/1ml
VDAC1 Polyclonal Antibody	Mitochondrial Protein Loading Control	ABP53121	30/100/200ul
Na ⁺ /K ⁺ -ATPase α1 Pab	Membrane Protein Loading Control	ABP51894	30/100/200ul
Anti-Plant Actin Mab (3T3)	Plants Protein Loading Control	A01050	50ul/200ul/1ml/10ml
Anti-Rubisco (Large Chain) Mab (9Y6)	Plants Protein Loading Control	A01110	50ul/200ul/1ml/10ml

*** Biotin, HRP, FITC, Cy3, AbFluor™ Conjugated Loading Controls are available. More products and bulk orders, please contact Abbkine or local distributors.

Epitope & Tag Antibodies

Product Name	Description	Cat. No.	Size
Anti-Flag Tag Mab (1B10)	WB, IF, IP	A02010	50ul/200ul/1ml/10ml
HRP Conjugated Anti-Flag Tag Mab (1B10)	HRP Conjugated	A02015	50ul/200ul/1ml
Anti-Flag Tag Mab (1B10), AbFluor™ 488	Green Dye Conjugated	A02010A488	100ul
Anti-Flag Tag Mab (1B10), AbFluor™ 594	Red Dye Conjugated	A02010A594	100ul
Anti-GFP Tag Mab (3D3)	WB, IP	A02020	50ul/200ul/1ml/10ml
Anti-GFP Tag Mab (3D3), HRP Conjugate	WB	A02020HRP	50ul/200ul/1ml
Anti-GST Tag Mab (2A8)	WB	A02030	50ul/200ul/1ml/10ml
Anti-HA Tag Mab (4F6)	WB, IF, IP	A02040	50ul/200ul/1ml/10ml
HRP Conjugated Anti-HA Tag Mab (4F6)	HRP Conjugated	A02045	50ul/200ul/1ml
Anti-HA Tag Mab (4F6), AbFluor™ 488	Green Dye Conjugated	A02040A488	100ul
Anti-HA Tag Mab (4F6), AbFluor™ 594	Red Dye Conjugated	A02040A594	100ul
Anti-His Tag Mab (5C3)	WB, IF, IP	A02050	50ul/200ul/1ml/10ml
Anti-Myc Tag Mab (2D5)	WB, IF, IP	A02060	50ul/200ul/1ml/10ml
Anti-MBP Tag Mab (9Y5)	WB	A02070	50ul/200ul/1ml/10ml
Anti-mCherry Tag Mab (9D3)	WB, IF	A02080	50ul/200ul/1ml/10ml
Anti-RFP Tag Mab (9D1)	WB	A02120	50ul/200ul/1ml/10ml
Anti-S Tag Mab (9T10)	WB	A02130	50ul/200ul/1ml/10ml
Anti-V5 Tag Mab (11D5)	WB, IF, IP	A02170	50ul/200ul/1ml/10ml
Anti-VSV-G-Tag Mab (14D2)	WB, IF, IP	A02180	50ul/200ul/1ml/10ml

*** Biotin, HRP, FITC, Cy3, AbFluor™ Conjugated Loading Controls are available. More products and bulk orders, please contact Abbkine or local distributors.

Agarose/Magnetic Beads conjugated Tag Antibodies

Product Name	Description	Cat. No.	Size
Anti-Flag Tag Mouse Mab (1B10), Agarose	Flag-Agarose	A02010AGB	100ul/400ul/2ml
Anti-Flag Tag Mouse Mab (1B10), Magnetic Beads	Flag-Magnetic Beads	A02010MGB	100ul/400ul/2ml
Anti-HA Tag Mouse Mab (4F6), Agarose	HA-Agarose	A02040AGB	100ul/400ul/2ml
Anti-HA Tag Mouse Mab (4F6), Magnetic Beads	HA-Magnetic Beads	A02040MGB	100ul/400ul/2ml
Anti-His Tag Mouse Mab (5C3), Agarose	His-Agarose	A02050AGB	100ul/400ul/2ml
Anti-His Tag Mouse Mab (5C3), Magnetic Beads	His-Magnetic Beads	A02050MGB	100ul/400ul/2ml
Anti-Myc Tag Mouse Mab (2D5), Agarose	Myc-Agarose	A02060AGB	100ul/400ul/2ml
Anti-Myc Tag Mouse Mab (2D5), Magnetic Beads	Myc-Magnetic Beads	A02060MGB	100ul/400ul/2ml
Anti-V5 Tag Mouse Mab (11D5), Agarose	V5-Agarose	A02170AGB	100ul/400ul/2ml
Anti-V5 Tag Mouse Mab (11D5), Magnetic Beads	V5-Magnetic Beads	A02170MGB	100ul/400ul/2ml

*** Agarose/Magnetic Beads conjugated Tag Antibodies can be directly used in IP or Co-IP. More products and bulk orders, please contact Abbkine or local distributors.

*** Magnetic Beads conjugated Tag Antibodies should be used with Magnetic separator.

IPKine™ Secondary Antibodies

Product Name	Description	Cat. No.	Size
IPKine™ HRP, Goat Anti-Mouse IgG LCS	light Chain Specific	A25012	100/500ul
IPKine™ HRP, Mouse Anti-Rabbit IgG LCS	light Chain Specific	A25022	100/500ul
IPKine™ HRP, Goat Anti-Mouse IgG HCS	Heavy Chain Specific	A25112	100/500ul
IPKine™ HRP, Goat Anti-Rabbit IgG HCS	Heavy Chain Specific	A25222	100/500ul

*** IPKine™ secondary antibody for WB after IP; LCS: Light Chain Specific, for eliminating heavy chain interference; HCS: Heavy Chain Specific, for eliminating light chain interference.

Protein Extraction & Quantification

Product Name	Description	Cat. No.	Size
Protein Quantification Kit (BCA Assay)	BCA Assay	KTD3001	500/2000/5000T
Protein Quantification Kit (Bradford Assay)	Bradford Assay	KTD3002	500/2000/5000T
Protease Inhibitor Cocktail (100X)	Protease Inhibitor	BMP1001	1/5*1ml
ExKine™ Total Protein Extraction Kit	Total Protein Extraction	KTP3006	50/200T
ExKine™ Nuclear and Cytoplasmic Protein Extraction Kit	Nuclear & Cytoplasmic Protein Extraction	KTP3001	50/200T
ExKine™ Nuclear Protein Extraction Kit	Nuclear Protein Extraction	KTP3002	50/200T
ExKine™ Cytoplasmic Protein Extraction Kit	Cytoplasmic Protein Extraction	KTP3003	50/200T
ExKine™ Total Membrane Protein Extraction Kit	Total Membrane Protein Extraction	KTP3004	50/200T
ExKine™ Membrane and Cytoplasmic Protein Extraction kit	Membrane & Cytoplasmic Protein Extraction	KTP3005	50/200T

*** More Protein or Organelle Extraction product, please contact Abbkine or local distributors.

ECL & PAGE Gel Staining & Protein Marker

Product Name	Description	Cat. No.	Size
SuperLumia ECL HRP Substrate Kit	Picogram Level	K22020	30/100/200/500ml
SuperLumia ECL Plus HRP Substrate Kit	Femtogram Level	K22030	30/100/200/500ml
Protein Gel Flash Staining Kit	Protein Gel Flash Staining	K21010	100/500ml
Colorcode Prestained Protein Marker (10-180 kDa)	Protein Marker	BMM3001	100/2*250/10*250ul
Colorcode Prestained Protein Marker (15-130 kDa)	Protein Marker	BMM3002	100/2*250/10*250ul

HRP & AP Conjugated Secondary Antibodies

Product Name	Description	Cat. No.	Size
HRP Conjugated Streptavidin	HRP- Streptavidin	A21000	100/500ul
HRP AffiniPure Goat Anti-Mouse IgG(H+L)	HRP Goat Anti-Mouse	A21010	100/500ul
HRP AffiniPure Goat Anti-Rabbit IgG(H+L)	HRP Goat Anti-Rabbit	A21020	100/500ul
HRP AffiniPure Goat Anti-Rat IgG(H+L)	HRP Goat Anti-Rat	A21040	100/500ul
HRP AffiniPure Donkey Anti-Sheep IgG(H+L)	HRP Donkey Anti-Sheep	A21060	100/500ul
HRP AffiniPure Goat Anti-Guinea Pig IgG(H+L)	HRP Goat Anti-Guinea Pig	A21070	100/500ul
HRP AffiniPure Goat Anti-Chicken IgY(H+L)	HRP Goat Anti-Chicken	A21080	100/500ul
HRP AffiniPure Goat Anti-Swine IgG(H+L)	HRP Goat Anti-Swine	A21090	100/500ul
AP AffiniPure Goat Anti-Mouse IgG(H+L)	AP Goat Anti-Mouse	A21110	100/500ul

*** HRP=Horseradish Peroxidase; AP=Alkaline Phosphatas; More [Enzyme Conjugated Secondary Antibodies](#), please contact Abbkine or local distributors.

IFKine™ Fluorescent Secondary Antibodies

Product Name	Description	Cat. No.	Size
IFKine™ Green Donkey Anti-Mouse IgG	Donkey Anti-Mouse (Green)	A24211	100/500ul
IFKine™ Green Donkey Anti-Rabbit IgG	Donkey Anti-Rabbit (Green)	A24221	100/500ul
IFKine™ Green Donkey Anti-Goat IgG	Donkey Anti-Goat (Green)	A24231	100/500ul
IFKine™ Orange Donkey Anti-Mouse IgG	Donkey Anti-Mouse (Orange)	A24311	100/500ul
IFKine™ Red Donkey Anti-Mouse IgG	Donkey Anti-Mouse (Red)	A24411	100/500ul
IFKine™ Red Donkey Anti-Rabbit IgG	Donkey Anti-Rabbit (Red)	A24421	100/500ul

*** Donkey host and other species of serum/IgG absorbed make IFKine™ ideal choice for fluorescence staining, especially [fluorescence multiple labelling](#).



Dylight™ Fluorescent Secondary Antibodies

Product Name	Color	Cat. No.	Size
DyLight 350 AffiniPure Goat Anti-Mouse IgG(H+L)	Purple	A23010	100/500ul
DyLight 350 AffiniPure Goat Anti-Rabbit IgG(H+L)		A23020	100/500ul
DyLight 405 AffiniPure Goat Anti-Mouse IgG(H+L)	Blue	A23110	100/500ul
DyLight 405 AffiniPure Goat Anti-Rabbit IgG(H+L)		A23120	100/500ul
DyLight 405 AffiniPure Rabbit Anti-Goat IgG(H+L)		A23130	100/500ul
DyLight 488 AffiniPure Goat Anti-Mouse IgG(H+L)	Green	A23210	100/500ul
DyLight 488 AffiniPure Goat Anti-Rabbit IgG(H+L)		A23220	100/500ul
DyLight 488 AffiniPure Rabbit Anti-Goat IgG(H+L)		A23230	100/500ul
DyLight 488 AffiniPure Goat Anti-Rat IgG(H+L)	Yellow	A23240	100/500ul
DyLight 549 AffiniPure Goat Anti-Mouse IgG(H+L)		A23310	100/500ul
DyLight 549 AffiniPure Goat Anti-Rabbit IgG (H+L)		A23320	100/500ul
DyLight 549 AffiniPure Rabbit Anti-Goat IgG (H+L)		A23330	100/500ul
DyLight 549 AffiniPure Goat Anti-Rat IgG (H+L)	Orange	A23340	100/500ul
DyLight 594 AffiniPure Goat Anti-Mouse IgG (H+L)		A23410	100/500ul
DyLight 594 AffiniPure Goat Anti-Rabbit IgG(H+L)		A23420	100/500ul
DyLight 594 AffiniPure Rabbit Anti-Goat IgG(H+L)		A23430	100/500ul
DyLight 594 AffiniPure Goat Anti-Rat IgG(H+L)	Red	A23440	100/500ul
DyLight 649 AffiniPure Goat Anti-Mouse IgG(H+L)		A23610	100/500ul
DyLight 649 AffiniPure Goat Anti-Rabbit IgG(H+L)		A23620	100/500ul
DyLight 649 AffiniPure Rabbit Anti-Goat IgG(H+L)		A23630	100/500ul
DyLight 649 AffiniPure Goat Anti-Rat IgG(H+L)	Odyssey Near-infrared	A23640	100/500ul
Dylight 680 AffiniPure Goat Anti-Mouse IgG(H+L)		A23710	100/500ul
Dylight 680 AffiniPure Goat Anti-Rabbit IgG(H+L)	Odyssey Far-infrared	A23720	100/500ul
Dylight 800 AffiniPure Goat Anti-Mouse IgG(H+L)		A23910	100/500ul
Dylight 800 AffiniPure Goat Anti-Rabbit IgG(H+L)	A23920	100/500ul	

*** More [Dylight™ Fluorescent Secondary Antibodies](#), please contact Abbkine or local distributors.

Blocking/Normal Serum

Product Name	Description	Cat. No.	Size
Normal Goat Serum	ELISA, ICC, IF, IHC	BMS0050	2/10/10*10ml
Normal Rabbit Serum	ELISA, ICC, IF, IHC	BMS0090	2/10/10*10ml
Normal Donkey Serum	ELISA, ICC, IF, IHC	BMS0140	2/10/10*10ml
Normal Human Serum	ELISA, ICC, IF, IHC	BMS0060	2/10/10*10ml
Normal Mouse Serum	ELISA, ICC, IF, IHC	BMS0070	2/10/10*10ml
Normal Rat Serum	ELISA, ICC, IF, IHC	BMS0100	2/10/10*10ml
Normal Horse Serum	ELISA, ICC, IF, IHC	BMS0150	2/10/10*10ml

*** [Blocking Serum](#) is strongly recommended as a blocking reagent to reduce background from non-specific, conserved sequence and/or Fc receptor binding. More Blocking Serum, please contact Abbkine or local distributors.

To get more information, Please contact Abbkine or visit <https://www.abbkine.com>.

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.

